Rat N-ERC/Mesothelin Assay Kit - IBL

96 Well

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96 Well

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1. Hino O, Kobayashi E, Nishizawa M, Kubo Y, Kobayashi T, Hirayama Y, Takai S, Kikuchi Y, Tsuchiya H, Orimoto K, et al. Renal carcinogenesis in the Eker rat. *J Cancer Res Clin Oncol*. 1995;121(9-10):602-5.
2. Yamashita Y, Yokoyama M, Kobayashi E, Takai S, Hino O. Mapping and determination of the cDNA sequence of the Erc gene preferentially expressed in renal cell carcinoma in the Tsc2 gene mutant (Eker) rat model. *Biochem Biophys Res Commun*. 2000 Aug 18;275(1):134-40.
3. Robinson BW, Creaney J, Lake R, Nowak A, Musk AW, de Klerk N, Winzell P, Hellstrom KE, Hellstrom I. Mesothelin-family proteins and diagnosis of mesothelioma. *Lancet*. 2003 Nov 15;362(9396):1612-6.
4. Scholler N, Fu N, Yang Y, Ye Z, Goodman GE, Hellström KE, Hellström I. Soluble member(s) of the mesothelin/megakaryocyte potentiating factor family are detectable in sera from patients with ovarian carcinoma. *Proc Natl Acad Sci U S A*. 1999 Sep 28;96(20):11531-6.
5. Robinson BW, Creaney J, Lake R, Nowak A, Musk AW, de Klerk N, Winzell P, Hellstrom KE, Hellstrom I. Soluble mesothelin-related protein--a blood test for mesothelioma. *Lung Cancer*. 2005 Jul;49 Suppl 1:S109-11.
6. Hassan R, Bera T, Pastan I. Mesothelin: a new target for immunotherapy. *Clin Cancer Res*. 2004 Jun 15;10(12 Pt 1):3937-42.
7. Onda M, Willingham M, Nagata S, Bera TK, Beers R, Ho M, Hassan R, Kreitman RJ, Pastan I. New monoclonal antibodies to mesothelin useful for immunohistochemistry, fluorescence-activated cell sorting, Western blotting, and ELISA. *Clin Cancer Res*. 2005 Aug 15;11(16):5840-6.
8. Hassan R, Remaley AT, Sampson ML, Zhang J, Cox DD, Pingpank J, Alexander R, Willingham M, Pastan I, Onda M. Detection and quantitation of serum mesothelin, a tumor marker for patients with mesothelioma and ovarian cancer. *Clin Cancer Res*. 2006 Jan 15;12(2):447-53.
9. Shiomi K, Miyamoto H, Segawa T, Hagiwara Y, Ota A, Maeda M, Takahashi K, Masuda K, Sakao Y, Hino O. Novel ELISA system for detection of N-ERC/mesothelin in the sera of mesothelioma patients. *Cancer Sci*. 2006 Sep;97(9):928-32.
10. Maeda M, Hino O. Molecular tumor markers for asbestos-related mesothelioma: serum diagnostic markers. *Pathol Int*. 2006 Nov;56(11):649-54.
11. Nakaishi M, Kajino K, Ikesue M, Hagiwara Y, Kuwahara M, Mitani H, Horikoshi-Sakuraba Y, Segawa T, Kon S, Maeda M, Wang T, Abe M, Yokoyama M, Hino O. Establishment of the enzyme-linked immunosorbent assay system to detect the amino terminal secretory form of rat Erc/Mesothelin. *Cancer Sci*. 2007 May;98(5):659-64.
12. Kuwahara M, Takeda M, Takeuchi Y, Kuwahara M, Harada T, Maita K. Transforming growth factor beta production by spontaneous malignant mesothelioma cell lines derived from Fisher 344 rats. *Virchows Arch*. 2001 May;438(5):492-7.
13. Tsuchiya H, Tsuchiya Y, Kobayashi T, Kikuchi Y, Hino O. Isolation of genes differentially expressed between the Yoshida sarcoma and long-survival Yoshida sarcoma variants: origin of Yoshida sarcoma revisited. *Jpn J Cancer Res*. 1994 Nov;85(11):1099-104.
14. Hino O, Shiomi K, Maeda M. Diagnostic biomarker of asbestos-related mesothelioma: Example of translational research. *Cancer Sci*. 2007 Aug;98(8):1147-51.
15. Hagiwara Y, Hamada Y, Kuwahara M, Maeda M, Segawa T, Ishikawa K, Hino O. Establishment of a novel specific ELISA system for rat N- and C-ERC/mesothelin. Rat ERC/mesothelin in the body fluids of mice bearing mesothelioma. *Cancer Sci*. 2008 Apr;99(4):666-70.

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