



## 検証試験 ( Total AGT vs Intact AGT )

作成日 : 2019年8月29日

検体 ( 健常人 )	希釈倍率 4,000倍 (ng/mL)			希釈倍率 10,000 倍 (ng/mL)		
	27412 (Total AGT)	27742 (Intact AGT)	Intact/Total (ratio)	27412 (Total AGT)	27742 (Intact AGT)	Intact/Total (ratio)
Serum 42	36,356	16,176	44.5%	39,280	15,760	40.1%
Serum 43	25,632	14,544	56.7%	26,230	14,080	53.7%
Serum 44	25,704	13,820	53.8%	28,130	14,430	51.3%
EDTA Plasma 42	28,612	14,168	49.5%	30,050	14,890	49.6%
EDTA Plasma 43	25,284	13,780	54.5%	26,790	14,830	55.4%
EDTA Plasma 44	21,240	11,892	56.0%	22,740	11,100	48.8%