


# Endothelin

- 血管収縮活性の高い 21mer と低い 38 mer を個々に検出できます。
- 世界初のエンドセリン (1 - 31) 測定キットを開発し幅広いエンドセリン研究が可能です。
- 心疾患、高血圧、脳血管障害、皮膚の炎症など幅広い研究分野で活用できます。



## ■ エンドセリン関連製品一覧

※ 受注生産品

F/P: ホルマリン固定・パラフィン包埋切片

索引	分類	製品番号		製品名	適用・測定範囲	容量	価格
Big Endothelin-1	Kit	27168	Human	Big Endothelin-1 Assay Kit - IBL	0.78-100 pg/mL	96 Well	¥98,000
	Kit	27167	Rat	Big Endothelin-1 Assay Kit - IBL	0.78-100 pg/mL	96 Well	¥95,000
	抗体	18203	Anti- Human	Big Endothelin-1 (22-38)	IHC (F/P)	50 ug	¥48,000
	抗体	18205	Anti- Rat	Big Endothelin-1 (22-38) Rabbit IgG Affinity Purify	IHC (F/P)	5 ug	¥15,000
Endothelin	抗体	18204	Anti-	Endothelin (15-21) (C) Rabbit IgG Affinity Purify	IHC (F/P)	50 ug	¥48,000
Endothelin-1	Kit	27165		Endothelin-1 Assay Kit - IBL	0.78-100 pg/mL	96 Well	¥98,000
	Kit	27177		Endothelin-1 (1-31) Assay Kit - IBL	1.56-200 pg/mL	96 Well	¥98,000
	抗体	18201	Anti-	Endothelin-1 Rabbit IgG Affinity Purify	IHC (F/P)	50 ug	¥48,000
Endothelin-2 (1-31)	Kit	※ 17178		Endothelin-2 (1-31) Assay Kit - IBL	3.91-500 pg/mL	96 Well	¥98,000
Endothelin-3	Kit	27169		Endothelin-3 Assay Kit - IBL	0.78-50 pg/mL	96 Well	¥98,000
	抗体	18202	Anti-	Endothelin-3 Rabbit IgG Affinity Purify	IHC (F/P)	50 ug	¥48,000
Endothelin-A Receptor	抗体	16201 ①	Anti- Human	Endothelin-A Receptor (A405) Rabbit IgG Affinity Purify	WB, IHC (F/P)	100 ug	¥58,000
Endothelin-B Receptor	抗体	16205 ②	Anti- Human	Endothelin-B Receptor (B27) Rabbit IgG Affinity Purify	WB	100 ug	¥58,000
	抗体	10253 ③	Anti- Human	Endothelin-B Receptor (8Z11) Mouse IgG MoAb	WB, IHC (F/P, Microwave 処理)	100 ug	¥58,000

### 参考文献:

1. Okishima N. *et al.* Specific sandwich-type enzyme immunoassays for smooth muscle constricting novel 31-amino acid endothelins. *Biochem Biophys Res Commun.* 1999; 256 (1), 1-5.
2. Tane N. *et al.* The effects of Endothelin-1 on Flap Necrosis. *Annals of Plastic Surgery.* 1995; 35(4), 389-395.
3. Unoki H. *et al.* Low-density lipoproteins modulate endothelial cells to secrete endothelin-1 in a polarized pattern: a study using a culture model system simulating arterial intima. *Cell and Tissue Research.* 1999; 295 (1), 89-99.
4. Hiraki H. *et al.* Regular immunohistochemical localization of endothelin-1 and endothelin-B receptor in normal, hyperplastic and eoplastic human adrenocortical cells. *Pathology International.* 1997; 47 (2-3), 117-125.
5. Muramatsu M. *et al.* Chronic hypoxia augments endothelin-B receptor-mediated vasodilation in isolated perfused rat lungs. *American Journal of Physiology.* 1999; 276 (2pt1), L358-L364.
6. Jun J. Wu *et al.* Inhibition of in vitro enteric neuronal development by endothelin-3: mediation by endothelin B receptors. *Development.* 1999; 126 (6), 1161- 1173.
7. Shimada K. *et al.* Modulation by endothelin-1 of lipopolysaccharide-induced cyclooxygenase 2 expression in mouse peritoneal macrophages. *European Journal of Pharmacology.* 1999; 376 (3), 285-292.
8. Uemura H. *et al.* Immunoreactive Endothelin-1 in the Neural Lobe of the Rat Pituitary Following Hemorrhage and Dehydration. *Endocrine Journal.* 1994; 41 (6), 685-691.
9. Yanagisawa H. *et al.* Altered expression of endothelin-1 and endothelial nitric oxide synthase in the juxtaglomerular apparatus of rats with HgCl<sub>2</sub>-induced acute renal failure. *TOXICOLOGY LETTERS.* 1998; 98 (3), 181-188.

