

IBL Code No.	製品名			適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KO463	Anti-Mouse	Trpm2	Rabbit PoAb	ELISA WB, ICC IHC	0.25	100 μL	25 μg	¥49,000
88KO464	Anti-Mouse	Trpm7	Rabbit PoAb	ELISA	0.25	100 μL	25 μg	¥49,000
88KM060	Anti-Rat	TRPM8 (CMR1)	Rabbit PoAb	ICC	0.25	100 μL	25 μg	¥49,000
88KM018	Anti-Rat	TRPV1 (VR-1)	Rabbit PoAb	WB, IHC	0.1	50 μL	5 μg	¥49,000
88KM019	Anti-Rat	TRPV2 (VRL-1)	Rabbit PoAb	IHC	0.1	50 μL	5 μg	¥49,000
88KM119	Anti-Mouse	TRPV4	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KX596	Anti-Human	TSHR	Mouse MoAb (Clone No.2E6)	ICC, FCM	0.25	200 μL	50 μg	¥55,000
88KG136	Anti-Human	TSHR/LGR3	Mouse MoAb (Clone No.3A6)	FCM	0.25	200 μL	50 μg	¥55,000
88KB555	Anti-Human	UPF3A	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KB557	Anti-Human	UQCRC1	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KB556	Anti-Human	UQCRC2	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KG403	Anti-Mouse	Usp11	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KH004	Anti-Carp	Vitellogenin	MoAb (Clone No.1G2)	ELISA WB	1	100 μL	100 μg	¥40,000
88KH006	Anti-Medaka	Vitellogenin	MoAb (Clone No.3C1)	ELISA WB	1	100 μL	100 μg	¥40,000
88KH005	Anti-Mummichog	Vitellogenin	MoAb (Clone No.3E11)	ELISA WB	1	100 μL	100 μg	¥40,000
88KH007	Anti-Red Sea Bream	Vitellogenin	MoAb (Clone No.5A4)	ELISA WB	1	100 μL	100 μg	¥40,000
88KG407	Anti-Mouse	Zmat4	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KG115	Anti-Human	β3-AR	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥39,000

【ELISA Kit】

分野	IBL Code No.	製品名		価格 (税抜)
環境関連	88KH003	コイ ビテロジェニン (Vg) ELISAキット	Carp Vitellogenin ELISA Kit	¥90,000
	88KH049	ウズラ ビテロジェニン (Vg) ELISAキット	Japanese Quail Vitellogenin ELISA Kit	¥75,000
	88KH075	メダカ ビテロジェニン (Vg) ELISAキット	Medaka Vitellogenin ELISA Kit	¥80,000
がん関連	88KK073	尿中ジアセチルスベルミン測定用ELISAキット	Urinary Diacetylspermine ELISA Kit	¥75,000
	88KK123	尿中ジアセチルスベルミン測定用ELISAキット	Urinary Diacetylspermidine ELISA Kit	¥100,000
	88KG573	可溶性CD147 (EMMPRIN) 測定用ELISAキット	Soluble CD147 (EMMPRIN) ELISA Kit	¥108,000
腎臓機能	88KK135	尿中クレアチニン測定用ELISAキット	Urinary Creatinine ELISA Kit	¥70,000

【その他の関連製品】

IBL Code No.	製品名	適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KH001-A	AGEs-BSA	標準タンパク	1	1000 μg	1000 μL	¥70,000

トランスジェニック社の

研究用試薬

抗体 & ELISA Kits

IBL Code No.	製品名		適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格	
88KR099	Anti-Human	17β HSD8	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KH043	Anti-	3-DG-Imidazolone	Mouse MoAb (Clone No.JNH 27)	IHC	0.25	200 μL	50 μg	¥55,000
88KE028	Anti-Rat	4F2 Heavy Chain (4F2hc:CD98)	Rabbit PoAb	WB, IHC	0.1	250 μL	25 μg	¥49,000
88KE020	Anti-Human	4F2 Heavy Chain (4F2hc:CD98)	Rabbit PoAb	IHC	0.5	500 μL	250 μg	¥49,000
88KB554	Anti-Human	ABCF1	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KB483	Anti-Human	ABCG2	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KR057	Anti-Mouse	Adducin	Rabbit PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KR081	Anti-Chicken	Adenosine deaminase2	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KG114	Anti-Human	AdipoR1	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥39,000
88KX577	Anti-Human	ADORA3	Mouse MoAb (Clone No.1C9)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB484	Anti-Human	ADRB2	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB485	Anti-Human	ADRB3	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KG132	Anti-	AGE-1	Mouse MoAb (Clone No.7C1)	ELISA, WB	0.25	40 μL	10 μg	¥55,000
88KG122	Anti-	AGE-3	Mouse MoAb (Clone No.9D8)	ELISA, WB	0.25	40 μL	10 μg	¥55,000
88KG133	Anti-	AGE-4	Mouse MoAb (Clone No.14B5)	ELISA, WB	0.25	40 μL	10 μg	¥55,000
88KH001	Anti-	AGEs	Mouse MoAb (Clone No.6D12)	ELISA, WB, IHC	0.25	40 μL	10 μg	¥55,000
88KH001-02	Anti-	AGEs	Mouse MoAb (Clone No.6D12) Peroxidase conjugated	ELISA, WB, IHC	0.1	200 μL	20 μg	¥70,000
88KH001-01	Anti-	AGEs	Mouse MoAb (Clone No. 6D12) Biotin conjugated	ELISA, IHC	0.25	40 μL	10 μg	¥70,000
88KR068	Anti-Human	Ajuba	Rabbit PoAb	IP, ICC	0.25	100 μL	25 μg	¥49,000
88KB544	Anti-Human	AKAP8	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KC598	Anti-Human	AKR7A3	Mouse MoAb (Clone No.2B8)	ELISA, WB IHC	0.25	200 μL	50 μg	¥55,000
88KR063	Anti-Rat	Aldehyde Oxidase	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KB511	Anti-Human	ANXA7	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KR078	Anti-Human	APH-1a	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KB538	Anti-Human	ARHGAP4	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB551	Anti-Human	ARNT	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KE027	Anti-Mouse	asc-type Amino Acid Transporter1 (Asc-1)	Rabbit PoAb	WB, IHC	0.1	250 μL	25 μg	¥49,000
88KB465	Anti-Human	ATP1B3	Mouse PoAb	WB	0.5	100 μL	50 μg	¥32,000
88KB542	Anti-Human	ATP2C1	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KR051	Anti-Human	Aurora-A	Rabbit PoAb	WB, ICC	0.25	100 μL	25 μg	¥49,000
88KM105	Anti-Mouse	Autotaxin /ENPP2	Rabbit PoAb	IHC	0.1	250 μL	25 μg	¥49,000
88KM106	Anti-Rat	Autotaxin /ENPP2	Rabbit PoAb	WB	0.1	250 μL	25 μg	¥49,000
88KR093	Anti-Mouse	BNAS2	Rabbit PoAb	WB, IP ICC	0.25	100 μL	25 μg	¥49,000
88KR086	Anti-Human	Brm	Rabbit PoAb	WB, IP	0.25	100 μL	25 μg	¥49,000
88KB486	Anti-Human	C5AR1	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KY041	Anti-Rat	Ca2+/Calmodulin-dependant Protein KinaseII δ1-δ4	Rabbit PoAb	WB, IHC	1	200 μL	200 μg	¥49,000
88KB512	Anti-Human	CACNB1	Mouse PoAb	WB	0.5	100 μL	50 μg	¥32,000
88KB513	Anti-Human	CACNB3	Mouse PoAb	WB	0.5	100 μL	50 μg	¥32,000

IBL Code No.	製品名			適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KO455	Anti-Mouse	CacnB4	Rabbit PoAb	ELISA, WB, IP, ICC	0.25	100 μg	25 μg	¥49,000
88KB514	Anti-Human	CACNG1	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KK054	Anti-Human	Calnexin Δ185-520	Mouse MoAb (Clone No.2E9)	WB	0.25	200 μg	50 μg	¥55,000
88KR084	Anti-Human	Calpain6	Rabbit PoAb	WB, IHC	0.25	200 μg	50 μg	¥49,000
88KR089	Anti-Human	CapG	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KR066	Anti-Mouse	Cation-Chloride-Cotransporter 9	Rabbit PoAb	IHC	0.25	100 μg	25 μg	¥49,000
88KB489	Anti-Human	CCBP2	Mouse PoAb	WB, FCM	1	50 μg	50 μg	¥32,000
88KM107	Anti-Zebrafish	Ccd1	Rabbit PoAb	WB, ICC	0.25	100 μg	25 μg	¥49,000
88KB488	Anti-Human	CCR5	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KB501	Anti-Human	CCRL2	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KS130	Anti-Human	CD147 (EMMPRIN)	Mouse MoAb (Clone No.2G2)	IP, FCM	0.25	200 μg	50 μg	¥55,000
88KS131	Anti-Human	CD298 (Na/K-ATPase b3)	Mouse MoAb (Clone No.5G10)	IP, IHC, FCM	0.25	200 μg	50 μg	¥55,000
88KN141	Anti-Human	CD3	Mouse MoAb (Clone No.NU-T3)	IHC, FCM	0.25	200 μg	50 μg	¥55,000
88KS017	Anti-Human	CD36	Mouse MoAb (Clone No.NF189)	FCM	0.25	200 μg	50 μg	¥30,000
88KO601	Anti-Human	CD44ICD (intracellular domain)	Rabbit PoAb	WB, IP, ICC, IHC	0.25	800 μg	200 μg	¥49,000
88KS125	Anti-Human	CD54 (ICAM-1)	Mouse MoAb (Clone No.YUK11)	IP, FCM	0.25	200 μg	50 μg	¥55,000
88KS127	Anti-Human	CD59 (HRF20)	Mouse MoAb (Clone No.YUK1)	IP, FCM	0.25	200 μg	50 μg	¥55,000
88KB467	Anti-Human	CD69	Mouse PoAb	WB, FCM	1	50 μg	50 μg	¥32,000
88KS128	Anti-Human	CD71 (TFRC)	Mouse MoAb (Clone No.YUK9)	IP, FCM	0.25	200 μg	50 μg	¥55,000
88KS124	Anti-Human	CD9 (MRP-1)	Mouse MoAb (Clone No.6D11)	WB, IP, IHC, FCM	0.25	200 μg	50 μg	¥55,000
88KB487	Anti-Human	CD97	Mouse PoAb	WB, FCM	1	50 μg	50 μg	¥32,000
88KS129	Anti-Human	CD98hc (4F2)	Mouse MoAb (Clone No.WK4)	IP, FCM	0.25	200 μg	50 μg	¥55,000
88KB468	Anti-Human	CEACAM5	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KH025	Anti-	CEL	Mouse MoAb (Clone No.KNH-30)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KH025-01	Anti-	CEL	Mouse MoAb (Clone No.KNH-30) Biotin conjugated	ELISA, IHC	0.25	200 μg	50 μg	¥70,000
88KH025-02	Anti-	CEL	Mouse MoAb (Clone No.KNH-30) Peroxidase conjugated	ELISA, IHC	0.25	200 μg	50 μg	¥70,000
88KB526	Anti-Human	CESK1	Mouse PoAb	WB	0.43	116 μg	50 μg	¥32,000
88KR072	Anti-Human	Chemokine-like factor 1	Rabbit PoAb	IHC	0.25	100 μg	25 μg	¥49,000
88KX591	Anti-Human	CHRM5	Mouse MoAb (Clone No.3A6)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KR103	Anti-Mouse	CIBZ	Rabbit PoAb	WB	0.25	400 μg	100 μg	¥49,000
88KB515	Anti-Human	CLIC2	Mouse PoAb	WB	1	50 μg	50 μg	¥32,000
88KB520	Anti-Human	CLIC3	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000
88KH011	Anti-	CML	Mouse MoAb (Clone No.CMS-10)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KH011-01	Anti-	CML	Mouse MoAb (Clone No.CMS-10) Biotin conjugated	ELISA, IHC	0.25	200 μg	50 μg	¥70,000
88KH011-02	Anti-	CML	Mouse MoAb (Clone No.CMS-10) Peroxidase conjugated	ELISA, IHC	0.25	200 μg	50 μg	¥70,000
88KH024	Anti-	CML	Mouse MoAb (Clone No.NF-1G)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KH024-01	Anti-	CML	Mouse MoAb (Clone No.NF-1G) Biotin conjugated	IHC	0.25	200 μg	50 μg	¥70,000
88KH024-02	Anti-	CML	Mouse MoAb (Clone No.NF-1G) Peroxidase conjugated	IHC	0.25	200 μg	50 μg	¥70,000
88KB469	Anti-Human	CSF2RA	Mouse PoAb	WB, FCM	0.32	152 μg	50 μg	¥32,000
88KR065	Anti-Rat	Cubilin	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KB503	Anti-Human	CXCR6	Mouse PoAb	WB	0.5	100 μg	50 μg	¥32,000
88KB547	Anti-Human	CYB5R3	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000

IBL Code No.	製品名			適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KB531	Anti-Human	RGS7	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KO454	Anti-Mouse	Rim1	Rabbit PoAb	ELISA, WB, IP	0.25	100 μg	25 μg	¥49,000
88KB497	Anti-Human	ROM1	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KB519	Anti-Human	RSC1A1	Mouse PoAb	WB, IHC, IF	0.5	100 μg	50 μg	¥32,000
88KX576	Anti-Human	RXFP3	Mouse MoAb (Clone No.2H10)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KR061	Anti-Human	S14G-Humanin	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KY008	Anti-Human	S19 Ribosomal Protein	Rabbit PoAb	WB	Lyoph.	Lyoph.	100 μg	¥45,000
88KR100	Anti-Mouse	Semaphorin 3B	Rabbit PoAb	WB, IP, IHC	0.25	400 μg	100 μg	¥49,000
88KO401	Anti-Mouse	Semaphorin 4A	Mouse MoAb (Clone No.HIAT-2)	ELISA, IP, FCM	0.25	100 μg	25 μg	¥55,000
88KO599	Anti-Mouse	Semaphorin 4B	Rat MoAb (Clone No.TK-2)	ELISA, IP, FCM	0.25	100 μg	25 μg	¥55,000
88KO453	Anti-Mouse	Semaphorin 4D (CD100)	Mouse MoAb (Clone No.SK-3)	ELISA, IP, FCM	0.25	100 μg	25 μg	¥55,000
88KO402	Anti-Mouse	Semaphorin 7A	Mouse MoAb (Clone No.SKK-7)	ELISA, IP, FCM	0.25	100 μg	25 μg	¥55,000
88KB536	Anti-Human	SH2D3A	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000
88KR055	Anti-Human	Siah1	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KR056	Anti-Human	Siah2	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KB549	Anti-Human	SLC16A1	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000
88KG404	Anti-Mouse	Slc25a19	Rabbit PoAb	IHC	0.25	100 μg	25 μg	¥49,000
88KB560	Anti-Human	SLC29A2	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000
88KB546	Anti-Human	SMAD6	Mouse PoAb	WB, FCM	0.45	109 μg	50 μg	¥32,000
88KO608	Anti-Mouse	Sox9	Rabbit PoAb	ELISA, WB, ICC, IHC	0.25	200 μg	50 μg	¥49,000
88KB481	Anti-Human	SPN	Mouse PoAb	WB, IF, FCM	0.5	100 μg	50 μg	¥32,000
88KB498	Anti-Human	SSR1	Mouse PoAb	WB, FCM	0.5	100 μg	50 μg	¥32,000
88KB499	Anti-Human	SSR2	Mouse PoAb	WB, FCM	1	50 μg	50 μg	¥32,000
88KX575	Anti-Human	SSTR2	Mouse MoAb (Clone No.1B7)	ELISA, IHC	0.25	200 μg	50 μg	¥55,000
88KM108	Anti-Rat	SulfFP1/sulf-1	Rabbit PoAb	WB	0.25	200 μg	50 μg	¥49,000
88KM110	Anti-Rat	SulfFP2/sulf-2	Rabbit PoAb (Antigen Affinity Purified)	WB, ICC	0.25	100 μg	25 μg	¥49,000
88KM109	Anti-Rat	SulfFP2/sulf-2	Rabbit PoAb (Protein G Purified)	WB, ICC	0.25	100 μg	25 μg	¥49,000
88KB479	Anti-Human	TACSTD1	Mouse PoAb	WB, IF	0.5	100 μg	50 μg	¥32,000
88KR076	Anti-Human	Tau Phospho-Specific (Se416)	Rabbit PoAb	WB	0.25	100 μg	25 μg	¥49,000
88KR111	Anti-Human	Tensin 2	Rabbit PoAb	WB, IHC	0.25	200 μg	50 μg	¥49,000
88KB541	Anti-Human	TIAM2	Mouse PoAb	WB, FCM	0.45	111 μg	50 μg	¥32,000
88KO578	Anti-Mouse	TP (Thromboxane A2 receptor)	Rabbit PoAb	ELISA, IF	0.25	100 μg	25 μg	¥49,000
88KB553	Anti-Human	TRIM45	Mouse PoAb	WB, IF	0.45	109 μg	50 μg	¥32,000
88KM120	Anti-Mouse	TRPA1	Rabbit PoAb	WB, IP	0.25	200 μg	50 μg	¥49,000
88KO457	Anti-Mouse	Trpc3	Rabbit PoAb	ELISA	0.25	100 μg	25 μg	¥49,000
88KO456	Anti-Mouse	Trpc3	Rabbit PoAb	ELISA, WB	0.25	100 μg	25 μg	¥49,000
88KO459	Anti-Mouse	Trpc5	Rabbit PoAb	ELISA	0.25	100 μg	25 μg	¥49,000
88KO458	Anti-Mouse	Trpc5	Rabbit PoAb	ELISA, WB	0.25	100 μg	25 μg	¥49,000
88KO460	Anti-Mouse	Trpc6	Rabbit PoAb	ELISA, ICC, IHC	0.25	100 μg	25 μg	¥49,000
88KO461	Anti-Mouse	Trpm1	Rabbit PoAb	ELISA, WB	0.25	100 μg	25 μg	¥49,000
88KO462	Anti-Mouse	Trpm2	Rabbit PoAb	ELISA, ICC	0.25	100 μg	25 μg	¥49,000

IBL Code No.	製品名			適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KB540	Anti-Human	GRB7	Mouse PoAb	WB, FCM	0.26	192 μg	50 μg	¥32,000
88KC597	Anti-Human	GRHPR	Mouse MoAb (Clone No.1H1)	ELISA, WB IHC	0.25	200 μL	50 μg	¥55,000
88KC595	Anti-Human	GRHPR	Mouse MoAb (Clone No.7G1)	ELISA, WB	0.25	200 μL	50 μg	¥55,000
88KB516	Anti-Human	GTF3C2	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KG408	Anti-Mouse	Haao	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB548	Anti-Human	HFE	Mouse PoAb	WB, FCM	0.43	114 μL	50 μg	¥32,000
88KB534	Anti-Human	HLA-A	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB532	Anti-Human	HLA-DOB	Mouse PoAb	WB	0.45	109 μL	50 μg	¥32,000
88KB535	Anti-Human	HLA-DPB1	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KB530	Anti-Human	HLA-DQB1	Mouse PoAb	WB	0.5	100 μL	50 μg	¥32,000
88KB558	Anti-Human	HLA-DRB1	Mouse PoAb	WB	0.5	100 μL	50 μg	¥32,000
88KB529	Anti-Human	HLA-DRB3	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB528	Anti-Human	HLA-E	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KR077	Anti-Mouse	Hixb9	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KR088	Anti-Mouse	HMGb3	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KX565	Anti-Human	HRH4	Mouse MoAb (Clone No.2D5)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KO607	Anti-Mouse	Hsd3b	Rabbit PoAb	ELISA, WB ICC, IHC	0.25	200 μL	50 μg	¥49,000
88KX574	Anti-Human	HTR5A	Mouse MoAb (Clone No.2D3)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KR050	Anti-Human	Humanin	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KB471	Anti-Human	IFNGR1	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB533	Anti-Human	IFNGR2	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB472	Anti-Human	IL1R1	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB482	Anti-Human	IL1R2	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB473	Anti-Human	IL2RB	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB474	Anti-Human	IL4R	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB475	Anti-Human	IL6R	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KB476	Anti-Human	IL7R	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KR101	Anti-Mouse	Indoleamine-pyrrole 2, 3 dioxygenase (IDO)	Rabbit PoAb	WB	0.5	300 μL	150 μg	¥49,000
88KR070	Anti-Mouse	Initiate Factor 3	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KR079	Anti-Mouse	Initiate Factor 5	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB522	Anti-Human	KCNIP2	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KB524	Anti-Human	KCTD15	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KB525	Anti-Human	KCTD18	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KB527	Anti-Human	KCTD4	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KR074	Anti-Rabbit	Keratin 12	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KR096	Anti-Mouse	KF-1	Rabbit PoAb	ICC	0.25	100 μL	25 μg	¥49,000
88KO603	Anti-Human	Klotho	Rat MoAb (Clone No.M2076)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO604	Anti-Human	Klotho	Rat MoAb (Clone No.KM2119)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO605	Anti-Mouse	Klotho	Rat MoAb (Clone No.Rink12)	IP	0.25	200 μL	50 μg	¥55,000
88KO606	Anti-Mouse	Klotho	Rat MoAb (Clone No.Rink14)	IHC, IF	0.25	200 μL	50 μg	¥55,000
88KG139	Anti-Human	KOR	Mouse MoAb (Clone No.5E3)	WB, IP ICC, FCM	0.25	100 μL	25 μg	¥55,000
88KB477	Anti-Human	LAMP2	Mouse PoAb	WB, IHC IF, FCM	0.49	102 μL	50 μg	¥32,000
88KG405	Anti-Mouse	Leprot1	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KG137	Anti-Human	LGR4	Mouse MoAb (Clone No.2H6)	WB, FCM	0.25	200 μL	50 μg	¥55,000
88KB478	Anti-Human	LRP1	Mouse PoAb	WB, IHC IF	1	50 μL	50 μg	¥32,000

IBL Code No.	製品名			適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格
88KE023	Anti-Human	L-type Amino Acid Transporter1 (LAT1)	Mouse MoAb (Clone No.4D9)	IHC	0.2	100 μL	20 μg	¥55,000
88KE026	Anti-Human	L-type Amino Acid Transporter1 (LAT1)	Rabbit PoAb	IHC	0.1	250 μL	25 μg	¥49,000
88KT117	Anti-Human	Macrophage (CD68)	Mouse MoAb (Clone No.PM-1K)	IHC	0.25	200 μL	50 μg	¥55,000
88KT022	Anti-Human	Macrophage Scavenger ReceptorA (MSR-A:CD204)	Mouse MoAb (Clone No. SRA-E5)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KT118	Anti-Human	Macrophage Scavenger ReceptorA (MSR-A:CD204)	Mouse MoAb (Clone No. SRA-C6)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KT013	Anti-Human	Macrophage Surface Antigen	Mouse MoAb (Clone No.AM-3K)	IHC	0.25	200 μL	50 μg	¥49,000
88KT014	Anti-Rat	Macrophage/Dendritic Cell	Mouse MoAb (Clone No.RM-4)	IHC	0.25	200 μL	50 μg	¥49,000
88KT015	Anti-Rabbit	Macrophage/Monocyte	Mouse MoAb (Clone No.RbM2)	IHC	0.25	200 μL	50 μg	¥49,000
88KJ091	Anti-Rat	Malignant Fibrous Histiocytoma (MFH) Cells	Mouse MoAb (Clone No.A3)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KT042	Anti-Rat	MARCKS	Rat MoAb (Clone No.MAR11/2)	WB	0.25	200 μL	50 μg	¥55,000
88KX446	Anti-Human	MAS1	Mouse MoAb (Clone No.15A10)	ELISA, WB	0.25	200 μL	50 μg	¥55,000
88KR083	Anti-Mouse	Matrix Gla Protein	Rabbit PoAb	IHC	0.25	200 μL	50 μg	¥49,000
88KR087	Anti-Human	MDC1 (NFB1/KIAA0170)	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KK053	Anti-Human	MEF	Rabbit PoAb	WB	0.1	100 μL	10 μg	¥49,000
88KA009	Anti-Rabbit	Metallothionein	Mouse MoAb (Clone No.1A12)	WB	1	100 μL	100 μg	¥45,000
88KA009-01	Anti-Rabbit	Metallothionein	Mouse MoAb (Clone No.1A12) Biotin conjugated	WB	1	100 μL	100 μg	¥70,000
88KR085	Anti-Mouse	MOK protein kinase	Rabbit PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KB480	Anti-Human	MPL	Mouse PoAb	WB, FCM	0.47	106 μL	50 μg	¥32,000
88KB466	Anti-Human	MS4A1	Mouse PoAb	WB, IHC FCM	1	50 μL	50 μg	¥32,000
88KR082	Anti-	Mu Phage gp44	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KR067	Anti-Human	Mucosa Associated Lymphoid Tissue 1 (MALT1)	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KG409	Anti-Mouse	Myipf	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KO588	Anti-Rat	NCX2	Rabbit PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KO580	Anti-Dog	NCX1	Mouse MoAb (Clone No.6H3-2b)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO584	Anti-Dog	NCX1	Rabbit PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KO581	Anti-Rat	NCX2	Mouse MoAb (Clone No.2S1-6)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO585	Anti-Rat	NCX2	Rabbit (NZW) PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KO586	Anti-Rat	NCX2	Rabbit (JW) PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KO582	Anti-Rat	NCX3	Mouse MoAb (Clone No.3U3)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO587	Anti-Rat	NCX3	Rabbit PoAb	WB, IHC	0.25	100 μL	25 μg	¥49,000
88KB521	Anti-Human	NEDD4L	Mouse PoAb	WB	0.45	109 μL	50 μg	¥32,000
88KH036	Anti-	Nitrotyrosine	Mouse MoAb (Clone No.2H1)		0.25	400 μL	100 μg	¥45,000
88KX564	Anti-Human	NPPFR1	Mouse MoAb (Clone No.2A10)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB494	Anti-Human	NPY5R	Mouse PoAb	WB, FCM	0.35	139 μL	50 μg	¥32,000
88KO610	Anti-Mouse	Nr5a1 (Ad4BP/SF-1)	Rat MoAb (Clone No.1B1F10)	ELISA, WB ICC, IHC, IF	0.25	200 μL	50 μg	¥55,000
88KB523	Anti-Human	NUDT9	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KB562	Anti-Human	NUP43	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KO140	Anti-Rat	Oazin (Antizyme inhibitor 1)	Mouse MoAb (Clone No.HI-12)	WB, ICC	0.25	200 μL	50 μg	¥55,000
88KE038	Anti-Human	Organic Anion Transporter1 (OAT1)	Rabbit PoAb	IHC	0.1	250 μL	25 μg	¥49,000
88KE031	Anti-Human	Organic Anion Transporter 2 (OAT2)	Rabbit PoAb	IHC	0.1	250 μL	25 μg	¥49,000

IBL Code No.	製品名		適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格	
88KE032	Anti-Human	Organic Anion Transporter 3 (OAT3)	Rabbit PoAb	IHC	0.1	250 μL	25 μg	¥49,000
88KE035	Anti-Rat	Organic Anion Transporter 3 (OAT3)	Rabbit PoAb	WB, IHC	0.1	250 μL	25 μg	¥49,000
88KE033	Anti-Human	Organic Anion Transporter 4 (OAT4)	Rabbit PoAb	IHC	0.1	200 μL	20 μg	¥49,000
88KR071	Anti-Rat	Organic Cation Transporter 3	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB509	Anti-Human	OXER1	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KX594	Anti-Human	OXGR1	Mouse MoAb (Clone No.2D4)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB504	Anti-Human	OXGR1	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KX592	Anti-Human	OXTR	Mouse MoAb (Clone No.1B8)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB517	Anti-Human	P2RX4	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB518	Anti-Human	P2RX5	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KB510	Anti-Human	P2RY8	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KO579	Anti-Rat	PAG (Phosphate-Activated lutaminase)	Mouse MoAb (Clone No.120)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KO583	Anti-Rat	panNCX	Mouse MoAb (Clone No.2S1-pan)	WB, IHC	0.25	200 μL	50 μg	¥55,000
88KB543	Anti-Human	PDCD8	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KR059	Anti-Mouse	PDX1	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KM037	Anti-Human	PEDF	Mouse MoAb (Clone No.14K)	WB	0.25	200 μL	50 μg	¥55,000
88KH012	Anti-	Pentosidine	Mouse MoAb (Clone No.PEN-12)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KH012-01	Anti-	Pentosidine	Mouse MoAb (Clone No.PEN-12) Biotin conjugated	ELISA, IHC	0.25	200 μL	50 μg	¥70,000
88KH012-02	Anti-	Pentosidine	Mouse MoAb (Clone No.PEN-12) Peroxidase conjug	ELISA, IHC	0.25	200 μL	50 μg	¥70,000
88KI044	Anti-Human	PERIOD1	Rabbit PoAb	WB	1	200 μL	200 μg	¥49,000
88KI047	Anti-Rat	PERIOD1	Rabbit PoAb	WB	1	200 μL	200 μg	¥49,000
88KI045	Anti-Human	PERIOD2	Rabbit PoAb	WB	1	200 μL	200 μg	¥49,000
88KI046	Anti-Human	PERIOD3	Rabbit PoAb	WB	1	200 μL	200 μg	¥49,000
88KI048	Anti-Rat	PERIOD3	Rabbit PoAb	WB	1	200 μL	200 μg	¥49,000
88KR097	Anti-Mouse	Phospholipid scramblase 2	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥49,000
88KM112	Anti-Rat	Phosphorylated TRPV1 (VR-1)	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KR090	Anti-Arabidopsis	Photorropin 2	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KR095	Anti-Arabidopsis	Phototropin 1	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KR064	Anti-Rat	Podocalyxin	Rabbit PoAb	WB, IF	0.25	100 μL	25 μg	¥49,000
88KR052	Anti-Mouse	POEM/nephronectin	Rabbit PoAb	IHC	0.25	200 μL	50 μg	¥49,000
88KG113	Anti-Human	PPAR γ	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥39,000
88KG406	Anti-Mouse	Prai4	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB496	Anti-Human	PTGER3	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KH010	Anti-	Pyrraline	Mouse MoAb (Clone No.H12)	ELISA, IHC	0.25	80 μL	20 μg	¥55,000
88KH010-01	Anti-	Pyrraline	Mouse MoAb (Clone No.H12) Biotin conjugated	ELISA, IHC	0.25	80 μL	20 μg	¥70,000
88KH010-02	Anti-	Pyrraline	Mouse MoAb (Clone No.H12) Peroxidase conjugated	ELISA, IHC	0.25	80 μL	20 μg	¥70,000
88KG134	Anti-Human	RAGE	Mouse MoAb (Clone No.1C5)	ELISA, WB, ICC, FCM	0.25	200 μL	50 μg	¥55,000
88KH039	Anti-Human	RAGE	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥49,000
88KE030	Anti-Rat	Related to b0, + Amino Acid Transporter (rBAT)	Rabbit PoAb	WB, IHC	0.1	250 μL	25 μg	¥49,000
88KB539	Anti-Human	RGS6	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000

IBL Code No.	製品名		適用	濃度 (mg/mL)	容量 (μg)	容量 (μL)	価格	
88KE021	Anti-Mouse	Cystin/Glutamic Acid Transporter (xCT)	Rabbit PoAb	WB, IHC	0.1	250 μL	25 μg	¥49,000
88KR080	Anti-Mouse	Dab1	Rabbit PoAb	IP	0.25	100 μL	25 μg	¥49,000
88KB493	Anti-Human	DARC	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KR104	Anti-Mouse	DC-STAMP (FIND)	Rabbit PoAb	ICC	0.25	100 μL	25 μg	¥49,000
88KG410	Anti-Mouse	Dcun1d4	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB545	Anti-Human	DDX11	Mouse PoAb	WB, IF	1	50 μL	50 μg	¥32,000
88KC570	Anti-Human	DDX39	Mouse MoAb (Clone No.2E4)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KG142	Anti-Human	DOR	Mouse MoAb (Clone No.1E7)	WB, IP, ICC, FCM	0.25	100 μL	25 μg	¥55,000
88KO602	Anti-	DYKDDDDK	Mouse MoAb (Clone No.2H8)	WB, IP, ICC, IF, FCM	0.25	200 μL	50 μg	¥55,000
88KO602-S	Anti-	DYKDDDDK	Mouse MoAb (Clone No.2H8)	WB, IP, ICC, IF, FCM	0.1	50 μL	5 μg	¥7,600
88KR062	Anti-Rat	Dynein	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KB537	Anti-Human	EDA2R	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB490	Anti-Human	EDG2	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB491	Anti-Human	EDG3	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB500	Anti-Human	EDG6	Mouse PoAb	WB	1	50 μL	50 μg	¥32,000
88KB550	Anti-Human	EFEMP2	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB552	Anti-Human	EIF4ENIF1	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KB559	Anti-Human	EME1	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KO609	Anti-Mouse	Emx2	Rabbit PoAb	ELISA, WB, ICC, IHC	0.25	200 μL	50 μg	¥49,000
88KR092	Anti-Alligator	Estrogen Receptor α (ER α)	Rabbit PoAb	WB	0.25	100 μL	25 μg	¥49,000
88KB470	Anti-Human	F3	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KO571	Anti-Mouse	Fc α/μ R	Mouse MoAb (Clone No.TX57)	IP, IF, FCM	0.25	200 μL	50 μg	¥55,000
88KO572	Anti-Mouse	Fc α/μ R	Mouse MoAb (Clone No.TX61)	IP, IF, FCM	0.25	200 μL	50 μg	¥55,000
88KX448	Anti-Human	FFAR2	Mouse MoAb (Clone No.16E4)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB492	Anti-Human	FPRL1	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KR102	Anti-Mouse	FRG	Rabbit PoAb	WB, IP	0.25	400 μL	100 μg	¥49,000
88KX593	Anti-Human	FZD8	Mouse MoAb (Clone No.1C3)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KR058	Anti-Mouse	G5PR	Rabbit PoAb	WB, IP	0.25	100 μL	25 μg	¥49,000
88KH040	Anti-Human	Galectin-3	Rabbit PoAb	WB	0.25	400 μL	100 μg	¥49,000
88KR069	Anti-Rat	Ghrelin	Rabbit PoAb	IHC	0.25	100 μL	25 μg	¥49,000
88KB563	Anti-Human	GOSR1	Mouse PoAb	WB, IF	0.42	119 μL	50 μg	¥32,000
88KB561	Anti-Human	GOSR2	Mouse PoAb	WB, IF	0.5	100 μL	50 μg	¥32,000
88KX589	Anti-Human	GPR108	Mouse MoAb (Clone No.2H9)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KG138	Anti-Human	GPR120	Mouse MoAb (Clone No.2B6)	ICC, FCM	0.25	100 μL	25 μg	¥55,000
88KB495	Anti-Human	GPR143	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KX447	Anti-Human	GPR149	Mouse MoAb (Clone No.1D2)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KB508	Anti-Human	GPR151	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KX590	Anti-Human	GPR162	Mouse MoAb (Clone No.1B9)	ELISA, IHC	0.25	200 μL	50 μg	¥55,000
88KG116	Anti-Mouse	GPR40	Mouse MoAb (Clone No.G16)	ICC, FCM	0.25	100 μL	25 μg	¥60,000
88KB502	Anti-Human	GPR56	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000
88KX443	Anti-Human	GPR65	Mouse MoAb (Clone No.11F1)	ELISA, WB	0.25	200 μL	50 μg	¥55,000
88KX444	Anti-Human	GPR68	Mouse MoAb (Clone o.15B8)	ELISA, WB	0.25	200 μL	50 μg	¥55,000
88KB505	Anti-Human	GPR77	Mouse PoAb	WB, FCM	0.5	100 μL	50 μg	¥32,000
88KB507	Anti-Human	GPRC5C	Mouse PoAb	WB, IF, FCM	1	50 μL	50 μg	¥32,000
88KB506	Anti-Human	GPRC5D	Mouse PoAb	WB, FCM	1	50 μL	50 μg	¥32,000