

Code No. 17316

Rat IgG

Volume : 1 mg

Lot No :

Introduction : Rat IgG is purified with protein A column from normal rat serum and it is lyophilized.
It does not contain preservative or stabilizer.

Purification : Affinity Purified with Protein A

Form : Lyophilized product in PBS

How to use : 1mL deionized water will be added to the product.

Stability : Lyophilized product, 5 years at 2 - 8°C
: Solution, 2years at -20°C

Application : This antibody can be used for IHC and WB as negative control.

For research use only, not for use in diagnostic procedures.