

Code No. 17317

Goat IgG

Volume : 1 mg

Introduction : The product is normal goat IgG purified from normal goat serum.

Purification : Purified with protein A column

Form : Lyophilized product in PBS (without additives)

How to use : 1.0 mL deionized water will be added to the product

Dilution : PBS (pH7.4)

Stability : Lyophilized product, 5 years at 2 - 8 °C
: Solution, 2 years at -20 °C

Application : Negative control in immunohistochemistry, western blotting or some tests.

For research use only, not for use in diagnostic procedures.