

Code No. 17508

**Anti-Mouse  
IgG (H+L) Rabbit IgG Fab'- HRP**

Volume : 0.2 mg

---

**Antigen** : Mouse IgG (H+L)

**Label** : HRP

**Purification** : Purified with antigen

**Form** : Lyophilized product from 1% BSA in PB (pH6.5)  
(containing ProClin<sup>®</sup>300 as preservative)

**How to use** : 1.0 mL deionized water will be added to the product

**Stability** : Lyophilized product, 5 years at -20 °C  
: Solution, 2 years at -20 °C

**Application** : This antibody can be used for western blotting in concentration of 10 - 100 ng/mL. (ECL)  
: This antibody can be used for ELISA or cell ELISA in concentration of 250 - 500 ng/mL.

**Specificity** : The antibody is absorbed by goat IgG, rabbit IgG and human IgG, but not absorbed by rat IgG.

---

*For Non-Clinical Research Use Only*