

Code No. 17562

**Anti-Rabbit  
IgG (H+L) Goat IgG Affinity Purify**

Volume : 1 mg

---

Antigen : Rabbit IgG (H+L)

Immunized animal :Goat

Purification : Affinity Purified with antigen

Process of absorption : This antibody was absorbed with Goat IgG, Mouse IgG, Rat IgG, Human IgG.

Specificity : Rabbit IgG : 100 %  
Mouse IgG · Rat IgG : less than 0.1 %  
Human IgG : less than 0.5 %

Form : Lyophilized product from PBS containing 1% BSA and 0.05% NaN<sub>3</sub>.

How to use : 1mL deionized water will be added to the product.

Application : This antibody can be used for IHC (1 - 10µg/mL) and WB (1 - 10µg/mL).

Method of Dilution : Dilute with PBS containing 1% BSA (pH 7.4)

Stability : Lyophilized product, 5 years at 2 - 8°C  
: Solution, 2years at -20°C