

Code No. 17567

**Anti-Human
IgG (Fc) Goat IgG Affinity Purify**

Volume : 1 mg
Lot No :

Antigen : Human IgG (Fc)

Immunized animal :Goat

Purification : Affinity Purified with antigen

Process of absorption : This antibody was absorbed with Goat IgG, Rabbit IgG, Mouse IgG, Rat IgG.

Form : Lyophilized product from PBS containing 1% BSA and 0.05% NaN₃.

How to use : 1mL deionized water will be added to the product.

Application : This antibody can be used for IHC (1 - 10μg/mL) and WB (1 - 10μg/mL).

Method of Dilution :Dilute with PBS containing 1% BSA (pH 7.4)

Stability : Lyophilized product, 5 years at 2 - 8°C
: Solution, 2years at -20°C

For Non-Clinical Research Use Only