

Code No. 17568

**Anti-Mouse
IgG (H+L) Rabbit IgG Affinity Purify**

Volume : 1 mg

- Antigen** : Mouse IgG (H+L)
- Purification** : Purified with antigen
- Form** : Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃
- How to use** : 1.0 mL deionized water will be added to the product
- Dilution** : PBS (pH7.4) containing 1% BSA
- Stability** : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C
- Application** : This antibody can be used for immunohistochemistry in concentration of 1 - 10 µg/mL.
: This antibody can be used for western blotting in concentration of 1 - 10 µg/mL
- Specificity** : The antibody is absorbed by goat IgG, rabbit IgG and human IgG, but is not absorbed by rat IgG.