

Code No. 17661

**Anti-Mouse  
IgG (H+L) Goat IgG Affinity Purify**

Volume : 1mg

---

**Antigen** : Mouse IgG (H+L)

**Purification** : Purified with antigen

**Form** : Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use** : 1.0 mL deionized water will be added to the product

**Stability** : Lyophilized product, 5 years at 2 - 8 °C  
: Solution, 2 years at -20 °C

**Application** : This antibody can be used for immunohistochemistry in concentration of 1 - 10 µg/mL.  
: This antibody can be used for western blotting in concentration of 1 - 10 µg/mL.

**Specificity** : The antibody is absorbed by rat IgG, goat IgG, rabbit IgG and human IgG.