

Code No. 18651

**Anti-Human
Thioredoxin (C) Rabbit IgG Affinity Purify**

Volume : 100 µg

Antigen	: Synthetic peptide in portion of C-terminus of Human Thioredoxin
Source	: Serum
Subclass	: IgG
Purification	: Purified with antigen peptide
Form	: Lyophilized product from 1% BSA in PBS containing 0.05% NaN ₃
How to use	: 1.0 ml distilled water will be added to the product
Dilution	: PBS (pH7.4) containing 1% BSA
Stability	: Lyophilized product, 5 years at 2 – 8 °C : Solution, 2 years at –20 °C
Application	: This antibody can be used for western blotting in concentration of 2 ~ 5 µg/ml.