

Code No. 50010

**Human CEA**

Volume : 0.1 mg

---

**Origin** : MKN-45 (Human gastric cancer cell line); supernatant

**Purity** : Purified by liquid chromatography  
: Immunological cross-reactive proteins such as NCA (non-specific cross reacting antigen) are removed by affinity purification.  
: Confirmed by western blotting

**Form** : Lyophilized product in PBS

**How to use** : Add 0.1 mL of deionized water to the product.

**Storage** : Lyophilized product, 2 - 8 °C  
Solution, -20 °C

---

*For Non-Clinical Research Use Only*

**Immuno-Biological Laboratories Co., Ltd.**  
1091-1 Naka Fujioka-Shi, Gunma 375-0005, JAPAN  
URL: <https://www.ibl-japan.co.jp/en/> E-mail: [do-ibl@ibl-japan.co.jp](mailto:do-ibl@ibl-japan.co.jp)