

Code No. 50010

Human CEA

Volume : 0.1 mg

Origin : MKN-45 (Human gastric cancer cell line); supernatant

Purity : Purified by liquid chromatography
: Immunological cross-reactive proteins such as NCA (non-specific cross reacting antigen) are removed by affinity purification.
: Confirmed by SDS-PAGE

Form : Lyophilized product in PBS

How to use : Add 0.1 mL of deionized water to the product.

Storage : Lyophilized product, 2 - 8 °C
Solution, -20 °C

For Non-Clinical Research Use Only