



Antibody Specificity Evaluation

#28065 Anti-Human GBF1 (T1337 Phosphorylated) Rabbit IgG Affinity Purify

#28063 Anti-Human GBF1 (1844) Rabbit IgG Affinity Purify



WB #28063 (Lot.0B-828)



W.B. #28065 (Lot: 0B-828)

| | |
|----------------|---------------------------|
| Date Reported: | 2 nd July 2020 |
| Conducted by | K. Ishikawa |
| Result: | Pass |

Test Details:

Cultured Huh7 cell line with stimulation 2-DG (-)/(+) and conducted W.B. using #28063 and #28065.

Comment:

The targeted bands were confirmed around expected size. As same as the reference (PMID: 18063581), 2-DG (+) lanes have stronger bands.

Lanes:

- 1: HuH-7 2DG(-)
- 2: HuH-7 2DG(+)
- 3: HuH-7 2DG(-) IP/#28063 (Lot: 0B-828)
- 4: HuH-7 2DG(+) IP/#28063(Lot: 0B-828)
- 5: HuH-7 2DG(-) IP/#28063 (Lot: 0B-806)
- 6: HuH-7 2DG(+) IP/#28063 (Lot: 0B-806)