

Freeze-Thaw Stability Test Result

[27900 Human Titin N-Fragment \(Urine\) ELISA Kit - IBL](#)

Sample Type: **Human Urine**

Sample		Result (pmol/L)	Dilution	Calc. (nmol/L)	%
No.1	FT0	1350.7	10	13.51	100
	FT1	1191.9	10	11.92	88.2
	FT2	1247.9	10	12.48	92.4
	FT3	1214.2	10	12.14	89.9
No.2	FT0	398.4	10	3.98	100
	FT1	356.1	10	3.56	89.4
	FT2	361.1	10	3.61	90.6
	FT3	371.9	10	3.72	93.3
No.3	FT0	293.8	10	2.94	100
	FT1	269.5	10	2.70	91.7
	FT2	274	10	2.74	93.3
	FT3	275.9	10	2.76	93.9

It is stable up to 3 times freeze-thaw repeat.

